

Efficient, Large-Scale Virus-Like Particle Manufacturing for Gene Editing by a GMP-Compliant Flow Electroporation® Platform

Isabel Daher¹, Tien Nguyen¹, Peter Gee¹, Sean Menarguez¹, Andrew Mancini¹, Marianna Romito¹, Ashley B. Strickland-Dietz¹, Lauren E. Unsworth¹, and James Brady¹
¹MaxCyte, Inc., Technical Applications, Rockville, MD

Abstract Virus-like particles (VLPs) are non-infectious, virally-derived nanoparticles that lack a viral genome. Upon packaging a CRISPR ribonucleoprotein complex (RNP), VLPs can target and transiently deliver cargo for genome editing *in vivo* to specific cell and tissue types. They are safer than traditional viral vectors as they are incapable of genomic integration. Scaling up VLP production can be cumbersome and time-consuming, requiring several rounds of optimization to transition from research to clinical scale. Additionally, consistency and reproducibility between batches is a major concern when relying on chemical methods of VLP production. Here, we utilized the MaxCyte ExPERT GTx®, a GMP-compliant electroporation instrument, to manufacture clinical-scale VLPs, in adherent and suspension HEK cells, packaged with CRISPR Cas9 or adenine base editor RNPs for genome editing in primary human cells. We found that electroporation consistently produced high yields of functional VLPs with an optimized electroporation workflow. We also demonstrated effective gene editing at several therapeutically relevant loci in primary hematopoietic cells, such as B2M and PD-1. Furthermore, the production of VLPs using electroporation exhibited favorable production kinetics compared to other transfection methods, enabling a one-day manufacturing process. Finally, we highlight the scalability of VLP production across a 400-fold volume range with minimal re-optimization, transfecting over one billion cells per production. In summary, our results show that MaxCyte's Flow Electroporation® technology is a viable means for consistent, efficient, and scalable manufacturing of VLPs for gene editing applications and has high promise to address the needs of future clinical and commercial manufacturing.

Figure 1. Experimental Workflow for VLP Production & Assessment of Yields




Figure 2. ExPERT™ Platform Enables Scalable VLP Production in Adherent HEK293 Cells




Figure 3. ExPERT™ Platform Enables Scalable VLP Production in Suspension HEK293 Cells




Figure 4. Electroporation Produces High Titer VLPs Sufficient for Primary Immune Cell Editing




Figure 6. Large-Scale VLP Production Requires Minimal Re-optimization on the ExPERT Flow Electroporation® Platform




Figure 5. Electroporation facilitates high incorporation of ABE8e into VLPs





Figure 7. ExPERT Platform Enables 400-Fold Volume Scalability



Summary

1. The MaxCyte ExPERT™ Flow Electroporation® platform enables VLP production, in both adherent and suspension HEK293 cells.
2. Incorporation of a variety of RNPs into VLPs via electroporation compares favorably to chemical transfection.
3. VLPs produced via electroporation are potent, directing high-efficiency gene editing at low doses, and can be achieved within a 1-day production timeline.
4. VLP production is scalable, from 5 million to 10 billion cells, without the need for additional process development, maintaining a titer comparable to small-scale.

Corresponding Author: James Brady; jamesb@maxcyte.com MaxCyte, Inc., Tel: (301) 944-1700 info@maxcyte.com • www.maxcyte.com

© 2025 MaxCyte, Inc. All Rights Reserved. MaxCyte®, ExPERT ATx®, ExPERT STx®, ExPERT GTx®, and Flow Electroporation® are trademarks of MaxCyte, Inc. registered in the U.S. Patent and Trademark Office. ExPERT®, ATx®, GTx®, and STx® are trademarks of MaxCyte, Inc. PATENT PENDING.

References: 1) An M, et al. Nat Biotechnol. 2024. doi: 10.1038/s41587-023-02078-y. PubMed PMID: 38191664. 2) Banskota S, et al. Cell. 2022;185(2):250-65.e16. doi: 10.1016/j.cell.2021.12.021. PubMed PMID: 35021064. 3) Hamilton JR, et al. Cell Rep. 2021;35(9):109207. doi: 10.1016/j.celrep.2021.109207. PubMed PMID: 34077734. 4) Kluesner MG, et al. Crispr J. 2018;1(3):239-50. doi: 10.1089/crispr.2018.0014. PubMed PMID: 31021262. 5) Raguram A, Banskota S, Liu DR. Cell. 2022;185(15):2806-27. doi: 10.1016/j.cell.2022.03.045. PubMed PMID: 35798006.