


General Incubator with Built-in Roller or Shaker

NB-205Q / NB-205QR


NB-205QR is the general incubator with built-in Roller Apparatus. NB-205Q is the incubator with built-in shaker.




Window
on Outer Door




Automatic Stop
when Opening Door



Brushless Magnetic
Induction Drive
(205Q)



Stainless Steel
& Round Corner



4ea of 2L
(205QR)



Multiple Indicators
& Alarm System


Specification

NB-205Q Only



MODEL		NB-205Q	NB-205QR
Chamber Volume		142L	
Built-in Apparatus		1 Shaker	2-layer Roller Rack (2 Bottles for Each Layer)
Temperature	Range	Ambient +5°C to 60°C	
	Accuracy	±0.25°C at 37°C	
	Increment	0.1°C	
	Control	Microprocessor Digital PID	
RPM	Range	30rpm to 300rpm	1rpm to 30rpm
	Accuracy	±0.1rpm	
	Increment	1rpm	0.1rpm
	Control	Microprocessor Digital Controller	
Heating		5-side Direct Heating	Forced Convection
Time	Range	Continuous or 99h 59min	
	Accuracy	0.1%	
	Increment	1 minute	
	Memory	-	Running time
Chamber Material		SUS304	DC Motor
Motor		Plate Type BLDC Motor	
Shelves		2 Shelves	1 Shelf
Door		Magnetic Door with Double-Pane Glass Window	
Display		LED Display	
Operating Panel		Tact Switch Button	
Dimensions	Interior	480(W)x435(D)x694(H)mm	
	Exterior	545(W)x545(D)x870(H)mm	
Weight		65kg	70kg
Power (110V available)		AC220V, 50/60Hz	AC220V, 50/60Hz
Power Consumption		410W	





Accessories



Spring Rack



Microplate Rack

100ml
x 16ea250ml
x 9ea500ml
x 5ea

Shelves