

Medium Shaker


NB-T101M / NB-T101MRC / NB-T101MRCS

NB-T101M




This is Compact Size Shaker with benefits such as small foot print, easy to move. Built-in Plate Type BLDC motor provides low noise, low vibration. Despite compact size, shaking is powerful and work place is large to load various vessels.

NB-T101MRC



NB-T101MRCS




Specification



MODEL		NB-T101M	NB-T101MRC	NB-T101MRCS
Motion		Orbital		
RPM	Range	30 to 300rpm		
	Accuracy	±1rpm		
	Increment	1rpm		
Orbital Diameter		22mm		
Time	Range	Continuous or 99h 59min 59s		
	Accuracy	±1%		
	Increment	1min		
Motor		Plate Type BLDC Motor		
Display		3.5-inch LCD Color Touchscreen		
Dimensions	Platform	460(W)x455(D)mm with Rubber Pad	460(W)x455(D)mm with Rubber Pad	420(W)x410(D)mm with Rubber Pad
	Exterior	465(W)x540(D)x125(H)mm	465(W)x540(D)x125(H)mm	454(W)x430(D)x127(H)mm
Weight		29kg	35kg	23kg
Power (110V available)		AC220V, 50/60Hz		AC220V, 50/60Hz
Power Consumption		20W		20W

Accessories



Spring Rack



Microplate
Rack x 6ea




2 Stairs



Tube Rack



100ml
x 16ea



250ml
x 9ea



500ml
x 5ea



Holders