

Roller Apparatus Max-Roll


NB-204MR / NB-204GR

NB-204MR, Remote Roller, is highly flexible in facilitating its operation not only incubator inside but also else where on work table.

NB-204GR is designed for using in the production of massive attachment cells. The rollers are available in 10 bottles per one deck and 100 bottles in 10 decks of roller.


Roller's Micro Stepping Motor makes less vibration and low noise for stable rotation.

NB-204MR



* NB-204MR can be used in NB-203XL or bigger model.


NB-204GR3



NB-204GR10



* NB-204GR can be used in NB-205QRXXXL / 203QRXXXL



Specification

MODEL		NB-204MR	NB-204GR3	NB-204GR10
Capacity (2L Bottle Length Below 250mm)		110mm to 120mm, 4ea	110mm to 120mm, 30ea	110mm to 120mm, 100ea
RPM	Range	1 to 30rpm (*Based on 2L Bottle RPM)		
	Accuracy	±1rpm	±0.1rpm	±0.1rpm
	Increment	0.1rpm	0.1rpm	0.1rpm
	Control	Microprocessor Digital Controller		
Time	Range	Continuous or 99h 59min		
	Accuracy	±1%		
Alarm	Function	-	Buzzer, Display	
Motor	Type	DC Motor		
	Control	32 bit DSP PWM Drive		
	Input Voltage	±10% DC24V		
Display		LED Display		
Operating Panel		Tact Switch Button		
Operating Environment	Temperature	0 to 50°C		
	Humidity	35 to 85%RH		
Dimensions	Roller Deck	362(W)x 299(D)mm	750(W) × 564(D)mm	750(W) × 564(D)mm
	Exterior	375(W)x 370(D)x345(H)mm	756(W)×625.5(D)×642(H)mm	756(W)×625.5(D)×1799(H)mm
Weight		12.5Kg	40kg	80kg
Power (110V available)		AC220, 50/60Hz	AC220, 50/60Hz	AC220, 50/60Hz
Power Consumption		13W	70W	75W